

POP-ON®



	SIZE		30/BOX REF	100/BOX REF
	Small	25mm	32301	32101
	Medium	29mm	32302	32102
	Intermediate	32mm	32303	32103
	Large	36mm	32304	32104
	X-Large	41mm	32305	32105

*Covered under Medicare
HCPCS code A4349*

TO SELECT THE PROPER CATHETER SIZE

Using the most appropriate catheter size will maximize comfort and wear time. For complete instructions on catheter sizing, print out the *Male External Catheter Users' Guide with Sizer* found under LITERATURE in this section of our web site. Be sure to print the sizer panel at 100% scale.

– OR –

Contact Rochester Medical to request a catheter sizer:
info@rocm.com

800-243-3315 USA TOLL-FREE
507-533-9309 PHONE