


ULTRAFLEX®

CE



SIZE		30/BOX REF	100/BOX REF
Small	25mm	33301	33101
Medium	29mm	33302	33102
Intermediate	32mm	33303	33103
Large	36mm	33304	33104
X-Large	41mm	33305	33105

*Covered under Medicare
HCPCS code A4349*

TO SELECT THE PROPER CATHETER SIZE

Using the most appropriate catheter size will maximize comfort and wear time. For complete instructions on catheter sizing, print out the *Male External Catheter Users' Guide with Sizer* found under LITERATURE in this section of our web site. Be sure to print the sizer panel at 100% scale.

– OR –

Contact Rochester Medical to request a catheter sizer:
info@rocm.com

800-243-3315 USA TOLL-FREE
507-533-9309 PHONE